



UTILITY PICK UP AND PRECISION DRESSING PLIERS

Easy To Use • Ergonomically Designed • Eco-Friendly

- Highly polished finish resists corrosion and increases instrument life throughout repeated sterilization
- High-quality hinge will not corrode or separate
- Serrated beaks for firm grip on materials during transport
- Guide pin assists in proper beak alignment
- Made of heavy gauge, heat-treated stainless steel to prevent tip bending and misalignment.

PRECISION DRESSING PLIERS

Used to grasp or transfer materials in and out of the oral cavity.

PART CODE:

DPU17

UTILITY PICK UP PLIERS

Used for transporting materials chairside and in the sterilization area.

PART CODES:

DP1	DP18
DP2	DP18L (locking)
DP3	DP19
DP17	

