



IMS[®] LILAC UTILITY GLOVES

Enhancing Safety In the Instrument Management Cycle

Lilac gloves are for use in the sterilization area and during disinfection processes in the operatory.

- Flocklined nitrile, lilac utility gloves do not contain latex and stand out to professionals as the color of dental hygiene
- Resist chemicals better than latex gloves
- Puncture resistant for improved protection against sharps injuries
- Molded texture improves tactile feel and grip
- Can be autoclaved up to 250°F/121°C up to five times



THE CLEANING PROCESS IS AN IMPORTANT STEP IN THE INSTRUMENT MANAGEMENT CYCLE. MAKE IT SAFER AND MORE EFFICIENT WITH HU-FRIEDY'S INSTRUMENT MANAGEMENT SYSTEM (IMS).

SIZE AND PART CODES

(three pairs per package)

Small	(Size 7)	40-060
Medium	(Size 8)	40-062
Large	(Size 9)	40-064
Extra Large	(Size 10)	40-066

OTHER IMS CLEANING & CARE PRODUCTS



ENZYMAX[®] ULTRASONIC CLEANING SOLUTIONS

- Dual-enzyme ultrasonic detergent specifically formulated for efficient ultrasonic cleaning or presoaking of cassettes and instruments
- Fast-acting, protects your instruments and eliminates the need for hand-scrubbing



IMS[®] DAILY CLEAN GENERAL PURPOSE DETERGENT

- Concentrated powder detergent for efficient ultrasonic cleaning or presoaking of cassettes and instruments
- Penetrating powder dissolves easily to efficiently remove blood, saliva and oral debris



SHINE reNEW[™] PRODUCTS

- Effectively remove stain and rust, restoring the luster to instruments and metals
- Quick and simple to use, eliminating the need for abrasive scrubs

Visit us online at Hu-Friedy.com

How the best perform

