




BLACK LINE RUBBER DAM CLAMPS

| An Exceptional Restoration Begins with Excellent Isolation

Hu-Friedy's 48 clamps, made in the USA from proprietary stainless steel, provide full-mouth coverage.

POINTS OF PERFORMANCE

- **Satin Steel®** which resists corrosion while aiding in flexibility and memory
- **A dark coating** allowing for improved intraoral contrast and visual acuity
- **Matte finish** for decreased light reflection and glare



ANTERIOR RUBBER DAM CLAMPS

00X
RDCM00X

Flat jaws and high bow for small lower canines and first premolars.

6X
RDCM6X

Slightly offset jaws for labial caries on anterior teeth.

9X
RDCM9X

Offset flat jaws and stiff bow for anterior teeth.

9SX
RDCM9SX


Offset flat serrated jaws for anterior teeth.

212X
RDCM212X

Flat jaws for labial caries on anterior teeth.

212SAX
RDCM212SAX

Slightly offset staggered jaws for odd shaped anteriors.



MOLAR RUBBER DAM CLAMPS (GENERAL)



8AX
RDCM8AX

Deeply festooned four-point jaws for molars.

8ADX
RDCM8ADX

Deeply festooned jaws and extended bow for distal access to molars.

WBAX
RDCMWBAX

Slightly festooned jaws for molars.

10X
RDCM10X

Flat jaws designed for small upper left and lower right molars; paired with RDCM11.

11X
RDCM11X

Flat jaws designed for small upper right and lower left molars; paired with RDCM10.

14X
RDCM14X

Deeply festooned jaws for partially erupted or irregularly shaped molars.

14AX
RDCM14AX

Deeply festooned jaws for partially erupted or irregularly shaped molars; larger than RDCM14.

W14AX
RDCMW14AX

Deeply festooned jaws for partially erupted or irregularly shaped molars.

18X
RDCM18X

Flat jaws and stiff bow for molars.

24X
RDCM24X

Flat jaws for upper right and lower left molars; paired with RDCM25.

25X
RDCM25X

Flat jaws for upper left and lower right molars; paired with RDCM24.

27X
RDCM27X

Flat jaws for large molars.

56X
RDCM56X

Slightly festooned jaws for large molars.

56SX
RDCM56SX

Slightly festooned serrated jaws for large partially erupted or irregularly shaped molars.

How the best perform




COMPLIMENTARY PRODUCTS



Rubber Dam Punch-Ainsworth
RDPA

Rubber Dam Forceps
RDF



Rubber Dam Frame, Child
RDCF5



Rubber Dam Frame, Adult
RDAF6

PREMOLAR RUBBER DAM CLAMPS

 0X RDCM0X Flat jaws and high bow for longer premolars.	 1X RDCM1X Slightly festooned jaws for upper premolars.				
 1AX RDCM1AX Deeply festooned jaws for canines and premolars.	 2X RDCM2X Flat jaws for premolars.			 2AX RDCM2AX Flat jaws for larger premolars.	 2ASX RDCM2ASX Flat serrated jaws for larger premolars.
 W2X RDCM2X Flat jaws for premolars.	 27NX RDCM27NX Slightly festooned jaws for premolars.			 29X RDCM29X Flat jaws for lower premolars.	 209X RDCM209X Small, flat jaws with stiff bow for lower molars.

UPPER MOLAR RUBBER DAM CLAMPS

	 4X RDCM4X Flat narrow jaws for small upper molars.	 W4X RDCMW4X Flat narrow jaws for small upper molars.	 5X RDCM5X Slightly festooned jaws for large upper molars.	 8X RDCM8X Deeply festooned jaws for upper molars.
	 26NX RDCM26NX Flat jaws for upper molars.	 30X RDCM30X Offset jaws for upper right molars; paired with RDCM31.	 31X RDCM31X Offset jaws for upper left molars; paired with RDCM30.	 201X RDCM201X Slightly festooned jaws for upper molars.

LOWER MOLAR RUBBER DAM CLAMPS

	 3X RDCM3X Flat narrow jaws for small lower molars.	 W3X RDCMW3X Flat narrow jaws for small lower molars.	 7X RDCM7X Flat jaws for lower molars.	 7AX RDCM7AX Slightly reversed festooned jaws for large lower molars.
	 W7X RDCMW7X Flat jaws for lower molars.	 12AX RDCM12AX Serrated offset jaws for lower right molars; paired with RDCM13A.	 13AX RDCM13AX Serrated offset jaws for lower left molars; paired with RDCM12A.	 26X RDCM26X Slightly reversed festooned jaws for small lower molars.

To learn more about Hu-Friedy's Black Line Rubber Dam Clamps, visit: www.Hu-Friedy.com/BlackLineRDC
 See all of Hu-Friedy's isolation and endodontic products at: www.Hu-Friedy.com/EndoEssentials
 Watch a short video highlighting the importance of isolation in dentistry by searching "Hu-Friedy Isolation" on YouTube

