

# Presidential® Forceps

Hu-Friedy Presidential® Forceps are superbly constructed from forged stainless steel, meticulously crafted and finished with a precision hinge.

## Lifetime Guarantee

All Hu-Friedy Presidential Forceps carry a lifetime warranty against defects in material and workmanship.


## Made in the U.S.A.

Hu-Friedy proudly handcrafts all Presidential Forceps in the U.S.A.

Serrations can be added to most Forceps upon request.



## Forceps Application



**FIRST BECAUSE WE LAST.**

## Forceps Application Reference Chart

	PRESIDENTIAL	INCISORS	CANINES	PREMOLARS	1st MOLARS	2nd MOLARS	3rd MOLARS	ROOTS	UPPER	LOWER
F1	●	●	●						●	
FAF1	●	●	●						●	
F10S	●				●	●	●		●	
F16	●				●	●				●
F16M	●				●	●				●
F17	●				●	●				●
F18L	●				●	●			●	
F18R	●				●	●			●	
F23	●				●	●				●
F53L	●				●	●			●	
F53R	●				●	●			●	
F62	●	●	●	●					●	●
F65	●	●						●	●	
F69	●							●	●	●
F88L	●				●	●			●	
F88R	●				●	●			●	
F99C	●	●	●	●					●	
F101	●			●						●
F150	●	●	●	●				●	●	
FAF150	●	●	●	●				●	●	
F150A	●	●	●	●					●	
F150AS	●	●	●	●					●	
F150S*	●	●	●	●	●	●		●	●	
F151	●	●	●	●				●	●	●
FAF151	●	●	●	●				●	●	●
F151A	●	●	●	●						●
F151AS	●	●	●	●						●
F151S	●	●	●	●	●	●		●	●	●
F203	●	●	●	●				●	●	●
F210H	●						●		●	
F210S	●						●		●	
F222	●						●			●
FAF222	●				●	●	●			●
F286	●	●		●				●	●	
F300	●							●	●	
F301	●							●	●	●
F150K*		●	●		●	●			●	
F151K*		●	●		●	●		●		●
F23K*					●	●				●
FAFX67A							●		●	
FAFX74N		●	●							●

Note: Forceps #62 & #101 may be used for deciduous or permanent teeth. \*Pedodontic.

# Hu-Friedy®

**FIRST BECAUSE WE LAST.**

www.hu-friedy.com • 1-800-HU-FRIEDY  
3232 N. Rockwell Street, Chicago, IL 60618

©2006 Hu-Friedy Mfg. Co., Inc. HU-FRIEDY and  
Presidential are registered trademarks of Hu-Friedy Mfg.  
Co., Inc., registered in the USPTO. All rights reserved.

HF-230/20M/03-06 R3

