


IMS 4-Step Process:

IMS is a self-contained system that simplifies the instrument management process. Your instruments are protected and secure so they last longer. Your staff benefits because instrument handling is decreased. Hu-Friedy's long lasting, low maintenance stainless steel cassettes hold each procedural set-up through the entire cleaning, sterilization, storage and presentation cycle. IMS was designed for you and your busy practice – it's adaptable, flexible, easy, safe and efficient.



1 CHAIRSIDE

Patient Preview & Presentation




2 CLEANING

Processed & Prepped



4 STORAGE

Organized & Ready for Use



3 STERILIZATION

Wrapped & Ready for Sterilization

UNIVERSAL ORTHODONTIC CASSETTE

Exterior Dimensions 7.25" x 5.75" x 3.25" (184mm x 146mm x 83mm)
 Holds 10 pliers, a Mathieu needle holder and 3 hand instruments in an accessory area.

Red	IMORTHO0
Gray	IMORTHO1
White	IMORTHO2
Orange	IMORTHO3
Yellow	IMORTHO5
Purple	IMORTHO6
Blue	IMORTHO8
Green	IMORTHO9



Hu-Friedy Mfg. Co., Inc.
 3232 N. Rockwell St.
 Chicago, IL 60618
 1-800-HU-FRIEDY
 www.hu-friedy.com

© 2005 Hu-Friedy Mfg. Co., Inc. Hu-Friedy is a registered trademark of Hu-Friedy Mfg. Co., Inc., registered in the USPTO. IMS, How Practices Make Perfect and Double Decker are trademarks/service marks of Hu-Friedy Mfg. Co., Inc. All rights reserved.

IMS-103/15M/11-05



Hu-Friedy® Orthodontic Cassettes



IMS™ Orthodontic Cassettes:

Hu-Friedy's Instrument Management System™ (IMS) is the most comprehensive, reliable and complete instrument processing and infection control system on the market today.

IMS combines sterilization, storage and organization of instruments in one integrated system. Use IMS to: improve staff and patient safety, maximize office efficiency, streamline instrument organization and protect cutters from dulling or chipping.

Electropolish:

Made of long-lasting electropolished stainless steel that provides a protective sheet against corrosion and a professional appearance for years to come

Thin Cassette Design:

Allows for more cassettes to be processed at one time, enhancing productivity

Durable Silicone Rails:


Designed to safeguard instruments from damage during processing

Efficient Hole Pattern:

Protects instruments from protrusion while allowing for enhanced cleaning and sterilization

Positive Locking Mechanism:

Provides a visual indication that cassette is locked for safety



Thin Orthodontic Cassettes:


- 40 - 50% thinner than standard cassettes allowing for more cassettes to fit into processing equipment.
- Ultra low silicone rails gently hold instruments in place during processing reducing the chance for damage.
- Available in two sizes.

SMALL THIN ORTHODONTIC CASSETTE

EXTERIOR DIMENSIONS 7" X 8" X 0.75" (178MM X 203MM X 19MM)
Holds up to 4 pliers, and 3 hand instruments.

IMORTHSM8

Blue



LARGE THIN ORTHODONTIC CASSETTE

EXTERIOR DIMENSIONS 8" X 11" X 0.75" (203MM X 280MM X 19MM)
Holds up to 7 pliers, and 4 hand instruments.

Blue

IMORTHLG8


Double-Decker Cassettes:

DOUBLE-DECKER™ SMALL ORTHODONTIC CASSETTE

Exterior Dimensions 4.5" x 8" x 1.5" (114mm x 203mm x 38mm)
Holds up to 4 pliers, a Mathieu needle holder and 2 hand instruments.

Red	IM9140-OR
Gray	IM9141-OR
White	IM9142-OR
Orange	IM9143-OR
Yellow	IM9145-OR
Purple	IM9146-OR
Blue	IM9148-OR
Green	IM9149-OR

*OPTIONAL INSTRUMENT RACK AVAILABLE:
Part Code IM914X-ORA



DOUBLE-DECKER™ LARGE ORTHODONTIC CASSETTE

Exterior Dimensions 6" x 8" x 1.5" (152mm x 203mm x 38mm)
Holds up to 6 pliers, a Mathieu needle holder and 3 hand instruments.

Red	IM5200-OR
Gray	IM5201-OR
White	IM5202-OR
Orange	IM5203-OR
Yellow	IM5205-OR
Purple	IM5206-OR
Blue	IM5208-OR
Green	IM5209-OR

*PICTURED WITH OPTIONAL INSTRUMENT RACK
Part Code IM520X-ORA

** ALSO AVAILABLE WITH AN ACCESSORY AREA
Part Code IM920X-OR

